



PT KIP MEDIKA INTERNASIONAL



PROTECTION SUIT

BETS DISPOSABLE

(Does not include Goggles, Shoe Covers, Masks and Gloves)

Number Medical Device Permit

KEMENKES RI AKD 21603121298



HEALTH EQUIPMENT PRODUCTION CERTIFICATE
LICENCE DEPARTMENT OF HEALTH RI
FK.01.02/1/1572-e/2021

HEALTH EQUIPMENT DISTRIBUTION CERTIFICATE
LICENCE DEPARTMENT OF HEALTH RI
FK.01.01/VI/4542-e/2021

DESCRIPTION :

This Personal Protective Equipment is specially designed for the needs of medical personnel who have direct contact with patients who have the communicable disease.

This Personal Protective Equipment is designed to provide protection for the wearer from water splashes, sprays or droplets that can penetrate.

COMPOSITION:

Microphorus
Water Resistant
Local / Domestic Products

STORAGE INSTRUCTION :

The product must always be tightly closed
Store in a room with good air ventilation
Keep away from sources of ignition,
incompatible substances and direct sunlight.

HOW TO USE :

1. Make sure the Hazmat Suit Jumpsuit Coverall is in good condition and not damaged
2. Wash hands with soap or use hand sanitizer in accordance with existing health standards
3. Wear shoes that can protect your feet
4. Then wear the Hazmat Suit Jumpsuit Coverall properly by inserting the legs first, attach the sleeves and close the Hazmat Suit Jumpsuit Coverall on the body by raising the zipper to the neck, leaving the hood or head protector of the Hazmat Suit Jumpsuit Coverall open behind the neck.



Produced and Distributed by:

PT. KIP MEDIKA INTERNASIONAL

Komplek Pergudangan Tritan Point Blok C-17

Jl. Raya Ketajen Desa Wedi, Kel. Wedi, Kec. Gedangan, Kab.

Sidoarjo, Prov. Jawa Timur - INDONESIA